**(DANISH)**  
**CINCH-90 INSTRUKTIONER**

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**(ENGLISH)**  
**CINCH-90 INSTRUCTIONS**

CINCH-90: a fast-setting, addition reaction silicone (polyvinyl siloxane) for use as an impression material in all crown and edentulous procedures. Available in automix cartridges only. Compatible with other polyvinyls. Contains no hazardous ingredients. May be disinfected in 2% glutaraldehyde solution.

Working Time @ 23°C (73°F)	60 seconds
Minimum Intra-oral Setting Time	90 seconds
Final Durometer	55 (Shore A)

**INSTRUCTIONS**

- Insert cartridge into impression gun, twist off sealing cap, express 5mm (1/4”) of material to insure equal flow. If apertures are clogged, clear blockage with an explorer before proceeding.
- Securely attach mixing tip to cartridge with 1/4 turn. Check for firm attachment before proceeding. Squeeze trigger with moderate, even pressure to extrude mix.
- Load tray. Inject additional material onto prepared teeth as “wash”.
- After extruding material, leave mixing tip attached until next use.

**IMPORTANT! Like other polyvinyls, Cinch-90 emits hydrogen after setting. For best results, delay pouring model for three (3) hours.**

**STORAGE AND USAGE**

- Store between 18°C and 24°C (64°F-75°F). Exposure to excess heat will shorten normal shelf-life (12 months).
- For best results, use at normal room temperatures. To extend work/set times, use slightly chilled material.
- Contamination by latex gloves can interfere with the setting reaction of polyvinyls. Polyvinyls are not compatible with polyether and polysulfide materials. Acrylic and methacrylate residues might inhibit set-ting of polyvinyls.

**ORDERING INFORMATION**

Stock No. S435 - Cinch Platinum Tubes, contains 1 tube base, 80ml, and 1 tube catalyst, 80ml.
Stock No. S454 - Cinch Platinunn Automix Cartridges, contains 2 - 50ml cartridges plus 2 mixing tips.
Stock No. S453 - Cinch-90 Automix Cartridges only, con-tains 2 - 50ml cartridges plus 2 mixing tips.
Stock No. S456 - Cinch Automomix Putty, cartridges only, contains 2 - 50ml cartridges plus 2 mixing tips.

**WARRANTY**  
Parkell win replace, free of charge, any material found to be defective for a period of 1 year from date of manufac-ture if handled and stored properly. Parkell does not accept liability for any loss or damage, direct or consequential, from the use of or the inability to use the product described herein. The user shall determine the suitability of the prod-uct’s intended use, and the user assumes all risk and liabil-ity in connection therewith. There is no implied warranty connected with the sale of this merchandise.

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