

INSTRUCTIONS FOR USE

DA: BRUGSANVISNING

DE: GEBRAUCHSANWEISUNG

EL: ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ

ES: MODO DE USO

FI: KÄYTTÖOHJE

FR: MODE D'EMPLOI

IT: ISTRUZIONI PER L'USO

NL: GEBRUIKSAANWIJZING

NO: BRUKSANVISNING

PT: INSTRUÇÕES DE UTILIZAÇÃO

SV: BRUKSANVISNING

Mfd. in Italy for:

parkell®

300 Executive Dr., Edgewood, NY 11717 USA 🏭
(800) 243-7446 • Parkell.com



Directa AB, Finvids väg 8,
SE-194 47 Upplands Väsby, Sweden

EN

EXPLANATION OF SYMBOLS USED

	Professional use only
	Follow instructions for use
	Temperature limit
	Keep away from sunlight
	Use-by date (expiration)

	Batch code
	Catalogue / stock number
	Package contents
	Do not use if package is damaged
	Manufacturer

IMPORTANT:

TriPhasix is a micro-fine particle powder alginate. It handles slightly different than traditional alginates (e.g. Jeltrate® etc.). By altering the amount of water (see below), you can mix a regular or light body to suit your needs. TriPhasix impressions are more like silicone. The surface is smooth, shiny and firm, with finer detail.

CHROMATIC FORMULA:

Violet/pink phase: Mix

Light Pink phase: Load tray

White phase: Seat Tray

High Compressive Strength: 1.2MPa

Flavor: French Vanilla

Setting time: 2 minutes, 40 seconds
from water being added.

DIRECTIONS FOR USE:

Immediately before taking the impression, have the patient rinse with water. Alginate works best when the surfaces are fully hydrated.

1. Apply alginate adhesive to tray (not included), or use a rim-lock tray.
2. For a quadrant impression, add one full measure of water to one full scoop of powder. Both measures should be filled to the top.
3. Spatulate till the mix reaches a smooth consistency and a homogeneous color. The color will change from violet/pink to light pink during the mixing process. At room temperature, Parkell alginate allows 45 seconds of mixing time. It sets in 2 minutes 40 seconds. (To speed setting, use warm water.)
4. Load the tray during the light pink phase and seat it in the mouth when it turns white. (If you want a heavy bodied material for extra drive, wait a few seconds after it turns white before seating the tray.)
5. After 2 minutes, 40 seconds from water being added, confirm set by touching material, and remove the impression.
6. Disinfect by spraying (use 1:10 dilution of sodium hypochlorite, 1:2 dilution of iodophor or 2% glutaraldehyde with phenolic buffer) Then rinse with warm water.
7. All alginates lose accuracy as they dry, so for greatest precision, pour as soon as possible. If not possible, seal in tub with wet paper towel and pour within 24 hours.

REFILLS:

Refills come with one canister and two sets of scoops for every 5lbs of alginate.

• Alumined moisture-proof envelope refills:

1 lb (REF S400), 5lbs (REF S401), 10lbs (REF S402), 25lbs (REF S404)

• Measuring Scoops for TriPhasix™ Alginate (REF S403)

WARRANTY AND TERMS OF USE:

For full Warranty and Terms of Use information, please see www.parkell.com.

Safety Data Sheets (SDS) are available at www.parkell.com.

Parkell's Quality System is certified to ISO 13485.