

Can I use the SEcure Primer as a bonding agent for other composites?

No. The Primer is specifically designed to work only with the SEcure resin cement. There is no reactivity between any other composite and the SEcure Primer.

How can I tell the difference between SEcure's 4:1 syringe mixing tip and my previous cement's 1:1 syringe mixing tip?

Since both the 4:1 and 1:1 tips have a brown hub it can be easy to mix them up. The 4:1 mixing tips that are required for successful dispensing of SEcure are slightly longer, and have an orange mixing mechanism within the clear component of the tip.

Which works better with the SEcure Resin Cement, the Self-Etching Primer, or Brush&Bond?

Brush&Bond and the SEcure Self-Etching Primer are equal in their ability to:
Increase the early bond strength between the tooth and the SEcure resin cement.
Desensitize and seal the prep.

Should the Primer be light cured prior to final seating of the crown?

No. The Primer should not be light cured.

Which component contains 4-META, the Self-Etching Primer, or the Resin Cement?

The 4-META is in the Self-Etching Primer.