

Rx Only •  S446






VINYL POLYSILOXANE TRAY ADHESIVE

**INSTRUCTIONS
FOR USE**

parkell®


DEVICE DESCRIPTION: To improve adhesion, wipe trays with alcohol before applying adhesive. You must allow to air dry completely before loading tray, about 5 minutes.

EXPLANATION OF SYMBOLS USED:

	Professional use only
	Follow instructions for use
	Catalogue / stock number
	Batch code
	Use-by date (expiration)

	Flammable Liquid
	Health Hazard
	Irritant (skin and eye)
	Aquatic Toxicity

A00218revB1119

Parkell, Inc., 300 Executive Dr., Edgewood, NY 11717 USA 
(800) 243-7446 • Parkell.com