Blu-Mousse®

VPS Bite Registration Material

A Legacy VPS Bite **Registration for Decades**

For decades, Blu-Mousse®, has been the bite registration material of choice for discerning dental clinicians worldwide. These astute dentists appreciate the ability of Blu-Mousse to precisely capture the interocclusal relationship, which enables labs to deliver crowns that often need minimal or no occlusal adjustment.

- Stays stacked on the occlusal surface during placement.
- ~ Optimal consistency allows patient to naturally close with zero resistance while capturing every detail.
- Hardens to a plaster-like final set with a durometer of 88 to ensure precise articulations.
- Does not distort between the time the bite is taken and when the lab articulates the models.
- Resists breakage on removal from the mouth, during trimming of ~ excess flash and when used to position models during articulation.

Set Times

✓ 30 Seconds | 60 Seconds | 2 Minutes

Additional Clinical Applications

- Scannable using latest technologies
- A matrix for fabricating direct and indirect provisionals.

SKU#	DESCRIPTION	SET TIME	SCENT	
S457S	Blu-Mousse Super-Fast	30 Sec.	Unscented	Blu-Mousse® and Green-Mousse® Kits include: • (2) 50ml cartridges (100ml ttl.) • (6) Mixing tips
S448S	Blu-Mousse Super-Fast	30 Sec.	Peppermint 🥭	
S440S	Blu-Mousse 60	60 Sec.	Unscented	
S438S	Blu-Mousse Classic	2 Min.	Unscented	
S437S	Blu-Mousse Classic	2 Min.	Peppermint 🥖	
S455S	Green-Mousse	2 Min.	Unscented	

ACCESSORIES & REFILLS

S445	Mixing Tips Refill — Pink Base (30 tips); for 50 ml Blu-Mousse and Green-Mousse Cartridges

S343 50ml Split Cartridge Dispensing Gun



Parkell





Also Available

Green-Mousse® VPS **Bite Registration Material**

- Greater flexibility
- 2-Minute set time
- 60-Durometer



