

# DIGITEST<sup>3</sup>

Digital Pulp Vitality Tester

## Hone in on the Source of Your Patient's Pain

- ✓ Single tooth testing
- ✓ Intuitive single button operation
- ✓ Multiple test settings to accommodate your personal preference
- ✓ Self-calibrating for quick easy setup
- ✓ Ergonomic design
- ✓ Four autoclavable probes (included) promote easy access to all tooth surfaces
- ✓ Standard 9V battery enables simple in-office replacement

## Elevate Your Diagnostic Skills to the Next Level


When percussion, ice, heat and radiographs are not enough, electric stimulation may be the diagnostic tool you need for determining the cause of the patient's discomfort. That's where the Digitest<sup>3</sup> comes in.



## How it Works

The Digitest 3 electrically stimulates the tooth to test for vitality. Comparing the response to the stimulus of the suspected tooth with a known healthy tooth's response to an identical stimulus can help you hone in on the source of the patient's pain.

With the ability to offer the additional diagnostic test of electric stimulation using the Digitest 3, dentists worldwide are routinely diagnosing and treating problematic teeth that otherwise might have to forgo treatment until a later date, when the problem more clearly presents itself.



SKU#	DESCRIPTION
------	-------------

D655	<b>Digitest 3 Pulp Vitality Tester;</b> (4) Autoclavable Probes; 1-standard long, 1-standard short, 1-labial, 1-lingual, (1) Ground Clip, (1) Lead Wire, (1) 9V Battery
------	---

ACCESSORIES	
-------------	--

D636	Precision Probe Set; (2) Probes; one labial, one lingual
D625	Standard Metal Probe Set; (2) Probes; one short & one long
D631GC	Pulp Tester Ground Clip & Lead Wire; (2) Clips, (1) Wire



For more information, visit [Parkell.com](http://Parkell.com) or contact your authorized Parkell distributor.

