

Features

TurboSensor+ 360
(SKU: D760)

TurboSensor
(SKU: D560-110)

| INSERT COMPATIBILITY | | |
|--|---|---|
| Any Parkell 30K Inserts | ✓ | ✓ |
| Other Market-Leading Brands of 30K Inserts | ✓ | |
| Any Parkell 25K Inserts | ✓ with Parkell 25/30K handpiece (sold separately) | ✓ |
| Other Market-Leading Brands of 25K Inserts | ✓ with Parkell 25/30K handpiece (sold separately) | |
| ENHANCED SERVICEABILITY | | |
| Plug-n-Play Accessories Detachable handpiece cord and footpedal assembly for quick in-office replacement, and reduced downtime when service needs happen. | ✓ | |
| CRUISE-CONTROL | | |
| Continuous scaling without foot fatigue Seamlessly initiated through smart-touch interface. | ✓ | |
| INFECTION CONTROL FEATURES | | |
| Detachable autoclavable handpiece | ✓ | |
| 1st-in-Class Biofilm Protection for Waterlines Built-in PureClarity™ waterline purification system for all scaler waterlines | ✓ | |
| Automated Purge Modes Preprogrammed to CDC guidelines (20 sec. and 2 min. auto cycles) | ✓ | |
| WARRANTY | | |
| 3-Year Warranty* <small>* Warranty for scaler unit itself. Different warranty periods apply to accessories to the scaler, including, but not limited to, detachable handpieces and footpedals.</small> | ✓ | ✓ |
| STYLE OPTIONS | | |
| Choose from a variety of colors to fit your style | ✓  |  |

rev1225



Learn more about the new TurboSensor+ 360 at [parkell.com](https://www.parkell.com)