



**SHOP  
HYGIENE**

Shop ultrasonic inserts, detachable handpieces, and more high-quality essentials by the trusted Parkell name at [parkell.com](http://parkell.com).

**PARKELL**

# Take Ultrasonic Scaling to the Next Level

*The new TurboSensor+ ultrasonic scaler meets the demands of today's dentists & hygienists*

Using ultrasonics for periodontal debridement saves time during hygiene appointments, increases productivity, and helps alleviate the hand fatigue often associated with the repetitive motions required by hand instrumentation. But it also presents some unique challenges, which means hygienists and dentists are constantly on the lookout for new technology that will not just simplify the preventive procedure, but offer additional protections for themselves and their patients. Parkell designed the new TurboSensor+ ultrasonic scaler to meet these evolving needs. In other words, the industry asked, and Parkell delivered.

"The TurboSensor+ is truly an exceptional scaler," said Parkell President David Mott. "It combines the quality and performance of its decades-trusted predecessor, TurboSensor, with features that have been highly requested by our loyal customers, most importantly, enhanced infection control and a truly open platform for flawless integration with other brand inserts."

**Remarkable Yet Familiar**

TurboSensor+ continues the tradition of being economically friendly to practices' bottom line thanks to its sound design and durability. While the new sensor maintains its predecessor's practical, low-profile layout and footprint, TurboSensor+ users will enjoy new and exciting features designed to enhance the scaling experience.

The heart of the TurboSensor+ is its IntelliSense processor, which gives clinicians the freedom to achieve consistent performance using any market-leading 30K insert. TurboSensor+ also simplifies sterilization protocols by offering a new detachable, autoclavable handpiece that's easy



to connect, disconnect, and reprocess. It's no secret that patients are becoming more aware of their safety while in the dental chair. With minimal additional effort, clinicians can ease their concerns by opening the sealed pouch and removing the autoclaved handpiece chairside.

**Scaling, Simplified**

TurboSensor+ has a superior water flow system for more precise and individualized aerosol management with every insert, as well as the ability to prime the autoclaved handpiece with fresh water following the removal of stagnant water during the purge process. Speaking of purging, TurboSensor+ includes foot pedal-controlled tap-on purging for simplified waterline cleansing, with 2 built-in timers for both the between-patient 30-second purge and the 2-minute purge at the start and end of the day. This eliminates worry over whether the waterline has been purged long enough.

TurboSensor+ offers easy setup thanks to a more compact and thought-out power cable assembly. A variety of vibrant colors are available, including pearl, tangerine, midnight blue, electric blue, or lilac for a finish that complements any clinical environment. Like all Parkell equipment, ultrasonic devices and inserts are designed, handcrafted, and serviced in the United States with a best-in-class 3-year warranty.

**FOR MORE  
INFORMATION:**

631.249.1134  
[parkell.com](http://parkell.com)