

AccuFilm® FASTCHECK

Pre-Mounted Articulating Film




Articulation. Simplified.

AccuFilm®, the world-renown brand recognized for accurately marking occlusion for decades introduces its simplest and thinnest film to date. At only 8 microns, AccuFilm® FastCheck incorporates the accuracy dental professionals routinely rely on with the added convenience of a built-in frame.

- ✓ Built-in frame facilitates easy positioning over the occlusal surfaces of posterior teeth without the need for forceps
- ✓ Incredibly thin film at only 8 microns consistently produces sharp, easy-to-interpret markings even on hard-to-test highly polished metal and ceramic surfaces
- ✓ Easily identifies the smallest occlusal discrepancies with precise marking on the points of contact
- ✓ Double-sided color coating for simultaneous marking of both arches
- ✓ Equally effective on both wet and dry metal, porcelain, and zirconia surfaces



Ultra-Thin!
Now only 8 microns



SKU#	DESCRIPTION
S054	AccuFilm® FastCheck Articulating Film Double-Sided— Blue/Blue; (100 strips)

See other side for our complete line of AccuFilm Articulating products →








For more information on AccuFilm FastCheck, visit Parkell.com, call 800-243-7446, or contact your authorized Parkell distributor.

parkell®




AccuFilm® Pre-Cut Film Canisters

SKU#	DESCRIPTION
 S014	AccuFilm® Marking Film Single-Sided – Black ; (1) canister of 280 pre-cut strips
 S015	AccuFilm® Marking Film Single-Sided – Red ; (1) canister of 280 pre-cut strips
 S017	AccuFilm® Marking Film Double-Sided – Red/Black ; (1) canister of 280 pre-cut strips
 S018	AccuFilm® Marking Film Double-Sided – Black/Black ; (1) canister of 280 pre-cut strips
 S021	AccuFilm® Marking Film Double-Sided – Red/Red ; (1) canister of 280 pre-cut strips


Single and double-sided film thickness is 19 microns.

AccuFilm® Booklet

SKU#	DESCRIPTION
 S053	AccuFilm® Marking Film Booklet Double-Sided – Red/Black ; (5) booklets (50 sheets/booklet, 250 sheets ttl.)

Double-sided film thickness is 19 microns.

AccuFilm® IV Brush-On Liquid

SKU#	DESCRIPTION
 S020	AccuFilm® IV Brush-On Marking Liquid; (1) 10ml bottle of marking liquid (Blue), (1) 20ml bottle of solvent

Double-sided film thickness is 19 microns.



For more information on AccuFilm articulating products, visit Parkell.com, call 800-243-7446, or contact your authorized Parkell distributor.

