

# Predicta® bulk

## BIOACTIVE

Dual-Cure Restorative



## Universal, Dual-Cure Composite for All Your Restorative Needs

A dual-cure composite that can be used to fill any size cavity preparation, all the way up to the occlusal surface, without any need for additional layering. From large Class V's to four surface posterior restorations, there are no limits to what you can do with Predicta® Bioactive Bulk.

- ✓ Predicta Bioactive Bulk is bioactive and releases calcium and phosphate ions and releases and recharges fluoride ions to stimulate mineral apatite formation and remineralization at the material-tooth interface
- ✓ Dual-cure for unlimited depth of cure
- ✓ Fluorescence for a more natural appearance
- ✓ Multiple viscosities and shades
- ✓ Highly radiopaque
- ✓ Indicated for Class I, II, III, IV, and V restorations
- ✓ Strong and durable to withstand occlusal forces
- ✓ Highly polishable
- ✓ Nanohybrid versatility






**A:** Large cervical lesions #'s 20 and 21 prepped  
**B:** Restored with Predicta Bulk

(Michael Karter, DDS)

## What is Bioactivity?

Known in the field to equate to stronger bonds between the restoration and tooth, penetration and filling of micro-gaps, reduction in sensitivity, guard against secondary caries, and the sealing of margins against microleakage.



	
<p><b>Regular Viscosity (High Viscosity)</b></p> <p>Ideally suited for clinical situations when placing larger restorations, where stability is desirable and material drift is not.</p>	<p><b>LV (Low Viscosity)</b></p> <p>Recommended for more conservative preps, where the ability to flow and self-level into tighter spaces is preferred for better adaptation of hard-to-access prep margins and line angles.</p>
<p>Both Regular and LV Predicta Bulk versions are appropriate for, and can be used in all restorative situations. These guidelines should be considered an aid when choosing a bulk-fill restorative that more closely matches your restorative needs and preferences.</p>	

SKU#	DESCRIPTION	
S615	A1/B1 Shade	<b>Predicta Bioactive Bulk-Fill Kits include:</b> <ul style="list-style-type: none"> <li>• (1) 5 ml hand-held syringe</li> <li>• (20) 17-gauge mixing tips</li> </ul>
S616	A2/B2 Shade	
S610	(LV) A1/B1 Shade	
S611	(LV) A2/B2 Shade	

### ACCESSORIES & REFILLS

S620	5 ml Dispensing Gun
S621	17-Gauge Dispensing Tip Refill (30 pcs)
S622	19-Gauge Dispensing Tip Refill (30 pcs)



For more information, visit [Parkell.com](http://Parkell.com) or contact your authorized Parkell distributor.

**parkell®**