





# Turbo SENSOR™

Ultrasonic Scaler



## The Backbone of Dental Hygiene

- ✓ Versatile and dependable
- ✓ Clearly marked intuitive controls
- ✓ Automatically adjusts frequency to power any 25KHz or 30KHz Parkell insert
- ✓ Autoclavable handpiece sheaths improve infection control
- ✓ On-demand foot-controlled Turbo Power Boost for heavy calculus
- ✓ Broad range of power easily accommodates both the sensitive patient with light deposits to the patient who hasn't had a professional scaling in years, as well as everyone in between
- ✓ Expanded low-power "Perio" mode

	Convenience and control with the TurboSensor's foot pedal
	Preinstalled in-line sediment water filter (REF D419). Unit comes with two total
	The male quick-connect coupler at the end of the water hose allows you to easily connect the unit to your water source.
	Includes two autoclavable handpiece sheaths, made of high-temperature resin to last and keep asepsis (REF D576).

“ Finding a piece of equipment as *durable, affordable, and effective* as the TurboSensor is like finding the Holy Grail of dentistry. ”

— Monica Rogers, BSDH, RDH

## Specifications:

- **Size of Unit:** 1½"H x 5½"W x 7½"D
- **Weight of Unit:** 2.70 lbs (1.22 kg) (43.2 oz)
- **Weight of Transformer:** 3.85 lbs (1.75 kg / 61.6 oz)
- **Length of Handpiece Cable:** 7' (215cm)
- **Length of Foot Pedal Cable:** 7' (215cm)
- **Length of Water Hose:** 7'6" (230cm)
- **Length of Power Cord:** 110V: 12'6" (380cm)
- **Total Length of Transformer Cable from Wall Outlet to Unit:** 12'6" (380cm)
- **Power:** 120VAC 50/60 Hz
- **Protection Against Electric Shock:** Class 1, Type B applied part
- **Protection Against Ingress of Liquids:**  
Foot Pedal & Scaler: IPX1 (drip proof); Power Supply: IPX0 (Ordinary)
- **Mode of Operation of Equipment:** Intermittent (10 minutes on, 5 minutes off)
- **Operating Conditions:** 15–30°C, 20–70% RH (non-condensing)
- **Certifications:** All electrical devices have been certified to IEC 60601-1, IEC 60601-1-2 and CAN/CSA C22.2 and are TUV listed
- 110V (220V/230V also available)

SKU#	DESCRIPTION
D560-110	<b>TurboSensor Ultrasonic Scaler Includes:</b> <ul style="list-style-type: none"> <li>• (1) Scaler</li> <li>• (2) Autoclavable Handpiece Sheaths</li> <li>• (1) In-Line Sediment Water Filter with Spare Disk</li> <li>• (1) Quick-Disconnect Water Connection</li> <li>• (1) Foot Controller</li> </ul>
<b>ACCESSORIES &amp; REFILLS</b>	
D419	Replacement Sediment Water Filters (10 disks)
D418	Water Filter Assembly Kit Includes: <ul style="list-style-type: none"> <li>• (2) Sediment Water Filters</li> <li>• (2) Clear Plastic Ferrules</li> <li>• (1) Female Connector</li> <li>• (1) Male Connector</li> </ul>
D576	Autoclavable Handpiece Sheath



For more information, visit [Parkell.com](http://Parkell.com) or contact your authorized Parkell distributor.

**parkell®**